

Bowl Stands

- Made of AISI 304 stainless steel
- Single and double bowl designs available
- Available with lower shelf
- Open frame design allows accessibility from any side
- Swivel castors
- Optional MR compatibility
- Four pcs 75 mm castors two with brake and two electrically conductive

Mechanical properties

Yield 0.2%-Rp 0.2 : 235
Yield 1%-Rp 1 : 265
Strength tensile-Rm : 550
Alongation A50 : 45

International Designations

EN 10088 Number : 1.4301
EN 10088 Name : X5CrNi18-10

Article Codes	L (mm)	W (mm)	H (mm)	Explanation	Weight(kg)
AB90.04	350	350	850	Single	6
AB90.05	350	700	850	Double	9



Technical Drawing

