

## Anaesthesia Screen, Fixed Size

- Model : OT60.03  
 Fixed width : Ø12x820x640 mm  
 Grossweight : 2 Kg  
 Packing Dimensions : 70x85x10cm
- Made of stainless steel
  - Including fixing clamp



## Anaesthesia Screen, Adjustable

- Model : OT60.04  
 Grossweight : 2 Kg  
 Packing Dimensions : 80x80x10 cm
- Made of stainless steel
  - Including fixing clamp
  - Adjustable height
  - Extendible to both sides



## IV Stand, 2 hooks

- Model : OT60.05  
 Grossweight Packing : 2 Kg  
 Dimensions : 15x130x15cm
- Made of stainless steel
  - Including fixing clamp
  - Adjustable height
  - 2 hooks



## IV Stand, 4 hooks

- Model : OT60.06  
 Grossweight Packing : 2 Kg  
 Dimensions : 35x125x35 cm
- Made of stainless steel
  - Including fixing clamp
  - Adjustable height
  - 4 hooks

