

## Shoulder Support

### General features:

Model	: OT60.09
Product explanation	: Shoulder support, pair
Supportmaterial	: Polyurethane
Material	: Stainless Steel(AISI 304, DIN X5 CrNi 18-10)
Grossweight	: 4 kg
Packing Dimensions	: 17x60x21cm
Product documentation	: CE
Classification	: Class I (MDD 93/42 EEC)

### Product Features

- Including fixing clamps
- Extendible
- Adjustable height
- Made of stainless steel



## Technical Drawing

