

Rectal Positioning Set

General features:

Model	: OT60.20
Product explanation	: Proctological Accessories
Material	: Stainless steel(AISI 304, DIN X5 CrNi 18-10)
Grossweight	: 11 kg
Packing Dimensions	: 61x76x26cm
Product Documentation	: CE
Classification	: Class I (MDD 93/42 EEC)

Product Features:

Body Material	: Ø16mm stainless rod
Support Material	: Polyurethane

- Set of rectoscopy roll and knee supports
- Including fixing clamps
- Rotatable and tiltable knee support
- Height and angle adjustable rectoscopy roll



Technical Drawing

