

Ophthalmology, Head and Wrist Support

General features:

Model	: OT60.21
Product explanation	: Ophthalmologic Head and Wrist Support
Material	: Stainless steel(AISI 304, DIN X5 CrNi 18-10)
Grossweight	: 9 kg
Packing dimensions	: 61x76x26cm

Body material	: Ø16mm Stainless Rod
Head support material	: Polyurethane

Product Features:

- Rotatable and tiltable head support and wrist support
- Height and angle adjustable
- Made of stainless steel
- Adjustable mounting from 400 mm to 470 mm



Angular movements of specific points

A : X - 360 degree
: Y - 60 degree
: Z - 60 degree

B : X - 360 degree
: Y - 60 degree
: Z - 60 degree

Technical Drawing

