

**Extension device****OT60.90****General features****Description**

- Orthopedic traction set is suitable for use with all INSPI TAL operating tables in orthopedic surgeries. Traction boots with up & down movement and rotation capability around its own axis.

<b>Model</b>	: OT60.90
<b>Product explanation</b>	: Orthopedic traction set
<b>Table connection features</b>	: Table body connection
<b>Material</b>	: Stainless steel(AISI 304, DIN X5 CrNi 18-10)
<b>Gross weight</b>	: 50 kg
<b>Packing dimensions</b>	: 65x150x50cm
<b>Product documentation</b>	: CE
<b>Classification</b>	: Class I (MDD 93/42 EEC)

**Product Features:**

<b>Body/Legmaterial</b>	: 30x50mm stainless box profile
<b>Supportmaterial</b>	: Ø18mm stainless rod
	: PU

**Angular movements of specific points**

<b>A</b>	: Y - 90 degree
<b>B</b>	: Y - 90 degree
<b>C</b>	: Y - 360 degree : Z - 350 degree
<b>D</b>	: X - 360 degree

This accessory is symmetric.

**Technical Drawing**